

GLOVE INTEGRITY TESTER

Glove Integrity Tester (multiple automatic)

For proper functionality and to ensure the air tightness of the Isolator, it is very important that the glove ports be inspected regularly and prior to use, for any possible leaks. The inspection of the glove ports is recommended to be performed with the help of the ***Glove Integrity Tester***.

The ***Glove Integrity Testers*** is an equipment for testing air tightness of the glove port system. It can test one (1) or four (4) gloves simultaneously, depending on the chosen model;

This equipment contains the followings components

- ◆ One mobile IsoCart:
 - Stainless steel
 - Stainless steel shelves
 - Casters
- ◆ One Pneumatic control panel:
 - Pressure regulators
 - Valves and adaptors
 - Mini Pressure gauge
 - Nylon flexible tubes and connectors
 - Mounting plate fabricated from stainless steel
- ◆ Four Seal plate per opening size:
 - Inflatable gasket fabricated from EPDM material
 - Gasket mounting plate fabricated from HDPE material
 - Differential pressure gauge type Dwyer Magnehelic
 - Nylon flexible tubes



4- Glove/Automatic Glove Tester