

# IsoTester-G

Multiple automatic  
Glove integrity tester

To ensure the air tightness and the proper functionality of any isolator

Ensures the integrity of up to four gloves and sleeves simultaneously.

HMI touchscreen

USB functionality to view results on a remote PC

Audit trail, regulations for electronic records and signatures

**The inspection of the gloves and sleeves is required to be done with the glove integrity tester**



PROTECTING LIFE  
PROTÉGER LA VIE  
PROTEGIENDO LA VIDA



## IsoTester - G BENEFITS

- Full testing cycle in under 5 minutes
- Ensures the integrity of both sleeves and gloves
- Can test individually, or up to four gloves simultaneously
- Creates positive pressure inside each glove and sleeve while monitoring the pressure via the PLC system
- Full Pdf report
- Easily identify leaks

### APPLICATIONS

THE GLOVE INTEGRITY TESTER IS AN EQUIPMENT FOR TESTING AIR TIGHTNESS OF A GLOVE PORT SYSTEM.

## Components

### Features

- > Siemens HMI Touchscreen and PLC
- > USB Data Logging, to view results on a remote PC
- > Can test up to 4 ports simultaneously
- > Full PDF report of system validation with all relevant data points and set times for each port
- > Pre-loaded recipes for differing glove/sleeve arrangements
- > Programmable custom parameters to suit specific needs
- > Audit trail, regulations for electronic records and signatures

### Mobile IsoCart

- > Stainless steel frame
- > Stainless steel shelves
- > Casters

### Requirements

- > Electrical: 115v/230v, 2A
- > Compressed air: 90-120psi

Testing plates can be customized for any glove port and any isolator (please inform us of the size of the glove ports)

### WARRANTY

- > 24 months

### DOCUMENTATION & SERVICES

- > Operation and maintenance manual
- > IQ/OQ sold separately