

PassBox

Clean room
transfer chamber

PURSAN9648US

Reduces transfer times
while also preventing
air convection between
different rooms

Minimize the pollution extent



Dynamic Pass Box



Passive Pass Box

PROTECTING LIFE
PROTÉGER LA VIE
PROTEGIENDO LA VIDA

PASSBOX

BENEFITS



- Reduces transfer times
- Prevents air convection between different rooms
- Minimize the pollution extent

APPLICATIONS

TRANSFER ITEMS TO DIFFERENT ROOMS WITH DIFFERENT CLEAN ROOM LEVELS EFFICIENTLY

TYPE	PASS BOX		AIR SHOWER PASS BOX	
MODEL	PB-01	PB-02	ASPB-01	ASPB-02
External Size (WxDxH) mm	660x570x580	760x670x680	660x570x1000	760x670x1100
Internal Size (WxDxH) mm	500x500x500	600x600x600	500x500x500	600x600x600
Material	304 stainless steel			
Interlock	Electronic			
UV Lamp	8W*1	15W*1	8W*1	15W*1
HEPA Filter	/	/	Efficiency ≥ 99.999% at 0.3μm	
Clean Air Speed	/	/	0.36 ~ 0.54m/s	
Gross Weight (kg)	55	70	80	100
Package Size (WxDxH) mm	780x690x770	880x790x870	780x690x1190	880x790x1290

DOCUMENTATION & SERVICES

- > Operation and maintenance manual

WARRANTY

- > 12 months