


IsoBooth™

Mobile modular
sampling & testing booth


Positive and negative pressure testing booth for COVID 19



Reversible glove ports for two sided operation

Safe & Flexible: With positive and negative pressure settings to safely test
from inside or outside the booth

PROTECTING LIFE
PROTÉGER LA VIE
PROTEGIENDO LA VIDA



ISOBOOTH™

BENEFITS

- Reversible glove ports for operation from inside or outside the booth.
- Positive & negative pressure to safely test patients with contamination control
- Protects health care workers and patients, offering a safer environment in crucial situations
- Uniform airflow across the filter face for efficient filtration
- High visibility windows allow for easy patient interaction
- No need to change PPE for the HCP after each test

APPLICATIONS

SAMPLING & TESTING OF INFECTIOUS PATIENTS, WITH CONTAINMENT CONTROL

AIR FILTRATION

- > **PRE-FILTER** Washable ASHRAE 30% wire frame
- > **FILTER** HEPA filtration (99.99% efficient at 0.3 µm)
- > **AIRFLOW** Uniform across filter face

DIMENSIONS (W x D x H)

> BOOTH DIMENSIONS

mm	1065 x 1264 x 2785 (Interior working height: 2174.5)
Inches	42" x 49.75" x 109.75" (Interior working height: 85.6")

ELECTRICAL CERTIFICATION

- > 120V

CONTROL PANEL

- > **FAN** On/Off - Speed control
- > **LIGHT** On/Off

MATERIALS OF CONSTRUCTION

- > **SIDE / TOP PANELS** Clear plastic 1/4" panels
- > **STRUCTURE** Lightweight aluminum extrusion
- > **DOOR** Clear panel

NOISE LEVEL

- > Quiet, less than 65 db

DOCUMENTATION & SERVICES

- > Operation and maintenance manual

WARRANTY

- > 12 months